Vaisala OPC UA Server Software

Setup and Configuration
for use with Vaisala viewLinc 5.1

EN

DE

FR

ES

PT

ZH

JA

VAISALA

PUBLISHED BY

Vaisala Oyj
Vanha Nurmijärventie 21, FI-01670 Vantaa, Finland
P.O. Box 26, FI-00421 Helsinki, Finland
+358 9 8949 1

Visit our Internet pages at www.vaisala.com.

© Vaisala 2020

No part of this document may be reproduced, published or publicly displayed in any form or by any means, electronic or mechanical (including photocopying), nor may its contents be modified, translated, adapted, sold or disclosed to a third party without prior written permission of the copyright holder. Translated documents and translated portions of multilingual documents are based on the original English versions. In ambiguous cases, the English versions are applicable, not the translations.

The contents of this document are subject to change without prior notice.

Local rules and regulations may vary and they shall take precedence over the information contained in this document. Vaisala makes no representations on this document's compliance with the local rules and regulations applicable at any given time, and hereby disclaims any and all responsibilities related thereto.

This document does not create any legally binding obligations for Vaisala towards customers or end users. All legally binding

obligations and agreements are included exclusively in the applicable supply contract or the General Conditions of Sale and General Conditions of Service of Vaisala.

This product contains software developed by Vaisala or third parties. Use of the software is governed by license terms and conditions included in the applicable supply contract or, in the absence of separate license terms and conditions, by the General License Conditions of Vaisala Group.

This product may contain open source software (OSS) components. In the event this product contains OSS components, then such OSS is governed by the terms and conditions of the applicable OSS licenses, and you are bound by the terms and conditions of such licenses in connection with your use and distribution of the OSS in this product. Applicable OSS licenses are included in the product itself or provided to you on any other applicable media, depending on each individual product and the product items delivered to you.


Table of contents

English.....	5
Deutsch.....	17
Français.....	31
Español.....	43
Português.....	57
中文.....	71
日本語.....	83

About Vaisala OPC UA Server

Vaisala OPC UA Server software supports the retrieval of real-time and historical data from Vaisala viewLinc Enterprise Server by third-party applications running on an OPC UA client machine.

When the Vaisala OPC UA Server receives a data request from an OPC UA client, it utilizes Vaisala API calls to retrieve data from viewLinc Enterprise Server. After authorizing the request, viewLinc Enterprise Server sends data via a secure transfer to the Vaisala OPC UA Server, which then transfers the data to the requesting OPC UA client.

System requirements

Availability	Dedicated server available 24 hours a day, 7 days a week
Server management	Recommended: Connected to an uninterruptible power supply (UPS)
	Recommended: Backup solution with support for open file backup
Operating System	Windows Server® 2019 (64-bit) Windows Server® 2016 (64-bit) Windows Server® 2012 R2 (64-bit) Windows® 10 Enterprise (64-bit)
Application disk space required	2 GB
Security certificate for web interface	Authorized TLS certificate and key
License key	OPC UA Server license key (on USB drive).

1) Vaisala OPC UA Server-signed certificate and key can be generated during installation.

License requirements

The Vaisala OPC UA Server software is shipped with a license key (found on the installation USB) to enable the OPC feature in viewLinc Enterprise Server. The matching license key must be entered in the Vaisala OPC UA Server Installation Wizard.

The OPC UA Server license must support an equal or greater number of devices as the viewLinc Enterprise Server license key.

To upgrade the Vaisala OPC UA Server software licensing size, you are only required to enter the new license key in viewLinc.

Security requirements

The Vaisala OPC UA Server can be set up with strict or relaxed levels of certificate authentication with viewLinc Enterprise Server, and strict or relaxed user authentication between Vaisala OPC UA Server and OPC UA clients. The decision to use strict or relaxed security depends on your company's security policy.

- **Certificate authentication:** If your security policy requires the use of trusted certificates, **strict** certificate authentication must be selected during installation. Certificate authentication refers to a security process which checks the level of security between viewLinc and the Vaisala OPC UA Server. **Relaxed** certificate authentication allows for the use of self-signed certificates.

For information on certificate requirements, see [OPC UA client certificate requirements \(page 7\)](#).

- **User authentication:** Strict user authentication makes sure that OPC UA client requests are only accepted from a recognized OPC UA client username/password. Relaxed user authentication automatically accepts anonymously generated requests from OPC UA clients.

The Vaisala OPC UA Server installation wizard provides the option to generate certificates for use with OPC UA clients.

Generate a Certificate	Install a Trusted Certificate
For companies with network access limited to a few PCs	For companies with remote network access needs
Created during Vaisala OPC UA Server installation	Certificate request generated internally and sent to a certificate signing authority for paid validation
Valid for up to 10 years	Valid for 2 years
Free	Cost varies/annual renewal fee
Use when relaxed security connections are allowed between Vaisala OPC UA Server and your OPC UA clients.	Use when strict security connections are required between Vaisala OPC UA Server and your OPC UA clients.

OPC UA client certificate requirements

OPC UA client certificates can be either self-signed certificates or CA-signed certificates. Both self-signed certificates and CA certificates must have the following minimum properties to be accepted.

1. **Key Usage:** the **keyUsage** property must have these options as a minimum:
 - **digitalSignature**
 - **nonRepudiation**
 - **keyEncipherment**
 - **dataEncipherment**
2. **Extended Key Usage:** the **extendedkeyUsage** property must have these properties:
 - **serverAuth**
 - **clientAuth**
3. **Subject Alternative:** the names list must have these entries as a minimum:
 - **URI.1** = Application URI
 - **DNS.1** = Fully qualified host name

Subject Alternative names list entry example:

- **URI.1** = **urn:Vaisala:OpcUaServer**
- **DNS.1** = **myhost.vaisala.com**

Allowing certificates to be trusted by the OPC UA Server

Both self-signed certificates and CA certificates must be made trusted in the OPC UA Server in order to be able to connect using them. For instructions on making the certificates trusted, see the following sections:

[Allowing self-signed certificates to be trusted \(page 7\)](#)

[Allowing CA-signed certificates to be trusted \(no certificate revocation list\) \(page 8\)](#)

[Allowing CA-signed certificates to be trusted \(with certificate revocation list\) \(page 8\)](#)

Allowing self-signed certificates to be trusted

If the OPC UA client has generated its own self-signed certificate, perform the following steps to allow the certificate to be trusted by the OPC UA server.

- ▶ 1. Allow the client to attempt a connection. Its certificate will be rejected, and the rejected certificate will be placed into the following folder:

```
<data location>\pki\DefaultApplicationGroup\rejected\certs
```

2. Move (cut/paste) the certificate from the **rejected** folder to the following location:

```
<data location>\pki\DefaultApplicationGroup\trusted\certs
```

3. Restart the OPC UA Server. The next connection attempt will be trusted.

Allowing CA-signed certificates to be trusted (no certificate revocation list)

If the OPC UA client certificate has been signed by a CA, perform the following steps to allow the certificate to be trusted by the OPC UA server.

These steps apply to CA certificates when a certificate revocation list is **not** associated with them. For instructions on allowing a CA certificate to be trusted when there is a certificate revocation list associated with it, see [Allowing CA-signed certificates to be trusted \(with certificate revocation list\) \(page 8\)](#).

- ▶ 1. Allow the client to attempt a connection. Its certificate will be rejected, and the rejected certificate will be placed into the following folder:

```
<data location>\pki\DefaultApplicationGroup\rejected\certs
```

2. Move (cut/paste) the certificate from the **rejected** folder to the following location:

```
<data location>\pki\DefaultApplicationGroup\trusted\certs
```

3. Copy the CA certificate to the following location:

```
<data location>>\pki\DefaultApplicationGroup\issuer\certs
```

The CA certificate must be in the **DER** format with a **.der** extension.

4. Make sure that **VOPCUAServer.cfg** has the following setting:
[opcua]
suppress_revoke_status_unknown=1
5. Restart the OPC UA Server. The next connection attempt will be trusted.

Allowing CA-signed certificates to be trusted (with certificate revocation list)

If the OPC UA client certificate has been signed by a CA, perform the following steps to allow the certificate to be trusted by the OPC UA server.

These steps apply to CA certificates when a certificate revocation list (CRL) is associated with them. For instructions on allowing a CA certificate to be trusted when there is **not** a CRL associated with it, see [Allowing CA-signed certificates to be trusted \(no certificate revocation list\) \(page 8\)](#).

- ▶ 1. Allow the client to attempt a connection. Its certificate will be rejected, and the rejected certificate will be placed into the following folder:

```
<data location>\pki\DefaultApplicationGroup\rejected\certs
```

2. Move (cut/paste) the certificate from the **rejected** folder to the following location:

```
<data location>\pki\DefaultApplicationGroup\trusted\certs
```

3. Copy the CA certificate to the following location:

```
<data location>>\pki\DefaultApplicationGroup\issuer\certs
```

The CA certificate must be in the **DER** format with a **.der** extension.

4. Copy the CRL to the following location:

```
<data location>\pki\DefaultApplicationGroup\issuer\crl
```

The CRL must be in the **DER** format with a **.crl** extension.

5. Make sure that **VOPCUAServer.cfg** has the following setting:
[opcua]
suppress_revoke_status_unknown=0
6. Restart the OPC UA Server. The next connection attempt will be trusted.

Installing Vaisala OPC UA Server

Before you start the Vaisala OPC UA Server Installation Wizard, make sure you have completed the installation prerequisite tasks outlined in the setup checklist.

Vaisala OPC UA Server software is installed as a Windows service. This service will only run on the LOCALSYSTEM user account, the Windows default account with full control of the system. If your computer system is configured with custom security settings, the Vaisala OPC UA services may need additional configuration to run. See [Troubleshooting \(page 14\)](#).

Table 1 Setup checklist

VOPC UA Server installation prerequisite tasks	
<input type="checkbox"/>	Vaisala OPC UA Server license key added to viewLinc. The same license key entered in viewLinc is entered during Vaisala OPC UA Server installation. Refer to the add license key instructions in the <i>viewLinc User Guide</i> .
<input type="checkbox"/>	A dedicated group and user configured in viewLinc. The dedicated group must have View permission to access the required Zone and/or Location data. The user account is used only to transfer data to the Vaisala OPC UA Server so that activity between viewLinc and the Vaisala OPC UA Server can be traced clearly in the Events log. Refer to the add group and add user instructions in the <i>viewLinc User Guide</i> .
<input type="checkbox"/>	viewLinc hostname identified, or IP address (a recognized certificate hostname is required if your network requires certificate authentication). ¹⁾
<input type="checkbox"/>	viewLinc port number identified (default is 443).
<input type="checkbox"/>	Vaisala OPC UA server port number identified (default is 55000).
<input type="checkbox"/>	Your OPC UA client username and password identified (if your network requires a unique user to authenticate data requests). ²⁾

- 1) *Your company security policy determines whether you are required to use strict or relaxed certificate authentication for connections between Vaisala OPC UA Server to viewLinc Enterprise Server*
- 2) *If strict user authentication is required to authorize data requests between Vaisala OPC UA Server and OPC UA clients, you must identify the username/password that will be used to initiate requests from the OPC UA client.*

Install Vaisala OPC UA Server software


Check your Group Security Policy settings to make sure your user account has the necessary permission to run software on the installation server.

- ▶ 1. On a dedicated server, insert the Vaisala OPC UA Server USB and run *VOPCUAServerSetup.exe*.

To install Vaisala OPC UA Server software on a remote server, copy the *VOPCUAServerSetup.exe* file from the USB to the destination server.

2. Select the installation language. This language setting is used in the wizard.
3. Make sure you have completed the setup prerequisites. Check each option to confirm, then click Next.
4. Enter your Vaisala OPC UA Server license key.
5. Accept the Vaisala General License Conditions.
6. Accept the default installation path for the software, or specify a new destination folder (location must have at least 2 GB of free disk space).
7. Accept the default installation path for data, or specify a new destination folder (location must have at least 1 GB of free disk space).
8. To configure the connection between viewLinc Enterprise Server and the Vaisala OPC UA Server, first choose the Certificate authentication settings:

Relaxed:

Data transfers permitted using a self-signed security certificate.

Strict:

Data transfers require that the viewLinc Enterprise Server provide a trusted certificate hostname. If this option is selected, the viewLinc hostname entered must match the trusted certificate, and the username and password must be configured in viewLinc.

9. Add the viewLinc connection details:

viewLinc IP address or hostname:

Type the viewLinc Enterprise Server hostname your network PCs use to connect with viewLinc through a browser, or the IP address. Note that a fully qualified domain name is required if you are setting up strict certificate authentication (for example, *viewLinc.mycompany.com*).

viewLinc Port number:

Accept the default port number, 443, or type a new port number.

viewLinc OPC UA username/password:

Type the dedicated username and password that was created in viewLinc. This user account is used by the Vaisala OPC UA Server to communicate with viewLinc.

10. Test the connection. When all viewLinc settings are valid, you can click Next to continue.

11. Choose to generate certificate files or upload trusted certificate and key files:

Keep existing certificate and key

For upgrade only. Choose this option to automatically use the certificate files currently installed on OPC UA clients.

Upload a certificate and key (trusted)

Choose this option if you already have trusted certificate and key files and they are available on your network.


Generate VOPC UA Server-signed certificate and key files

Choose this option if your security policy permits a relaxed security connection with OPC UA clients, or you want to purchase a trusted certificate at a later time.

You can use the generated **.csr** file to purchase a trusted certificate for installation on OPC UA clients. See [Updating certificates \(page 14\)](#).

12. To configure the connection between Vaisala OPC UA Server and your OPC UA client, choose the user authentication requirement:

Relaxed:

Ensures data transfers are accepted without requiring a dedicated username/password (anonymous data requests from an OPC UA client are permitted).

Strict

Data transfers are only recognized if the request is from a recognized user account. After selecting this option, enter the accepted username/password combination (see step 14).

13. Add the Vaisala OPC UA Server port number the OPC UA client will use to communicate with the Vaisala OPC UA Server (default is 55000).
14. If you selected strict user authentication, enter the designated OPC UA client username and password authorized to generate requests.
15. Click **Install**
16. Click **Finish** to complete the installation wizard.
17. To verify the installation was successful, open the Services window (open the Windows Start menu and type **Ser v i c e s**), and locate Vaisala OPC UA Server in the Name column. The status should read, **Running**. If the service is not running, see [Troubleshooting \(page 14\)](#).

viewLinc data parameters

After installation of the Vaisala OPC UA Server is complete, OPC UA clients can be configured to retrieve the following viewLinc parameters:

Parameter
Device name
Probe serial number
Device serial number
Device channel number
Calibration date
Location name
Location time stamp
Location value (measurement): realtime
Location value (measurement): historical
Location units

Maintenance activities

CAUTION! Changes to the VOPC UA Server configuration file should only be performed by system administrators.

Updating certificates

You can update an expired security certificate, or update the Vaisala OPC UA Server system to use a trusted certificate.

- ▶ 1. Stop the Vaisala OPC UA Server service.
2. Locate the certificate and key folders in the default data installation path, for example: `... \Public Documents \Vaisala \Vaisala OPC UA Server \pki \DefaultApplicationGroup \OWN \certs \ ... \Public Documents \Vaisala \Vaisala OPC UA Server \pki \DefaultApplicationGroup \OWN \private \`
3. Replace the existing `\certs` folder file, `application_rsa_sha256.der`, and `\private` folder file, `application_rsa_sha256_key.pem`.
4. Start Vaisala OPC UA Server service.

Changing security levels

Depending on your security policy, you may need to change security certificate or user authentication requirements. A full description of strict vs. relaxed authentication is covered in [Security requirements \(page 6\)](#).

- ▶ 1. Stop the Vaisala OPC UA Server service.
2. Start the installation wizard, `VOPCUAServerSetup.exe`.
3. To change the certificate authentication security level, go to **viewLinc Enterprise Server Connection** step and change the certificate authentication to **Strict** or **Relaxed**.
4. To change the user authentication security level, go to **Your OPC UA Client Connection** step and change the user authentication to **Strict** or **Relaxed**.
5. Complete the installation wizard steps as required.
6. Start the Vaisala OPC UA Server service.

Troubleshooting

Vaisala OPC UA Server Service not running

By default Vaisala OPC UA Server software installs as a Windows service which is permitted to run only on the LOCALSYSTEM user account; this is the Windows default account which has full control on the system. If your computer system is configured with custom security settings, the Vaisala OPC UA Server service may need additional configuration to run.

1. Open the Windows **Start** menu and open or type **Services**.

2. Locate the **Vaisala OPC UA Server Service** and open the **Properties** window.
3. On the **Log on** tab select **This account**, then enter an authorized username and password. The user must have read and execute permission on all installed program files sub folders, inheritable Full Control access on all sub folders, and network access rights to allow it to connect to other systems.

OPC UA Service fails to load

If the OPC UA Service fails to load, check the *Log\VOPCUAServerLog* file and the *Log\VOPCUAServer_trace.log* file for the following messages:

- Vaisala OPC UA Server failed to load Configuration. [error code]
- Vaisala OPC UA Server failed to create add policy for level 1. [error code]
- Vaisala OPC UA Server failed to create add policy for level 2. [error code]
- Vaisala OPC UA Server initialize failed. [error code]
- Vaisala OPC UA Server failed to create certificate store. [error code]
- Vaisala OPC UA Server failed to create certificate. [error code]
- Vaisala OPC UA Server failed to create certificate request. [error code]

The [error code] shown at the end of the message is an internal error code for support purposes. Contact Vaisala Support to investigate why initialization is failing.

OPC UA client not connecting

If you generated security certificate files during installation, you may want to install the generated **.crt** file on connecting OPC UA client machines. This prevents possible connection errors when not using a trusted certificate.

1. On each OPC UA client, copy the generated certificate file (*VOPCUAServer-CA.crt*) to any desktop location, then right-click on the file to select **Install Certificate**.
2. In the **Certificate Import Wizard Welcome** screen, select **Local Machine**.
3. On the **Certificate Store** screen select **Place all**, click **Browse**, and then select **Trusted Root Certification Authorities**. If you receive an unknown publisher warning, click **OK**.
4. Click **Finish**, then **Yes**.

Certificate rejected when attempting to connect

OPC UA client certificates can be either self-signed certificates or CA-signed certificates. Both self-signed certificates and CA certificates must have certain minimum properties in place, and they must be allowed to be trusted by the OPC UA server.

For information on correct certificate properties, see [OPC UA client certificate requirements \(page 7\)](#).

For instructions on allowing self-signed and CA-signed certificates to be trusted by OPC UA server, see the following sections:

- [Allowing self-signed certificates to be trusted \(page 7\)](#)
- [Allowing CA-signed certificates to be trusted \(no certificate revocation list\) \(page 8\)](#)
- [Allowing CA-signed certificates to be trusted \(with certificate revocation list\) \(page 8\)](#)

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[Allowing self-signed certificates to be trusted \(page 33\)](#)

[Allowing CA-signed certificates to be trusted \(no certificate revocation list\) \(page 34\)](#)

[Allowing CA-signed certificates to be trusted \(with certificate revocation list\) \(page 34\)](#)

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1. Allow the client to attempt a connection. Its certificate will be rejected, and the rejected certificate will be placed into the following folder:

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