Setup Guide

System Setup and Configuration

Vaisala viewLinc Monitoring System

VAISALA

PUBLISHED BY

Vaisala Oyj
Vanha Nurmijärventie 21, FI-01670 Vantaa, Finland
P.O. Box 26, FI-00421 Helsinki, Finland
+358 9 8949 1

Visit our Internet pages at www.vaisala.com.

© Vaisala 2019

No part of this document may be reproduced, published or publicly displayed in any form or by any means, electronic or mechanical (including photocopying), nor may its contents be modified, translated, adapted, sold or disclosed to a third party without prior written permission of the copyright holder. Translated documents and translated portions of multilingual documents are based on the original English versions. In ambiguous cases, the English versions are applicable, not the translations.

The contents of this document are subject to change without prior notice.

Local rules and regulations may vary and they shall take precedence over the information contained in this document. Vaisala makes no representations on this document's compliance with the local rules and regulations applicable at any given time, and hereby disclaims any and all responsibilities related thereto.

This document does not create any legally binding obligations for Vaisala towards customers or end users. All legally binding

obligations and agreements are included exclusively in the applicable supply contract or the General Conditions of Sale and General Conditions of Service of Vaisala.


This product contains software developed by Vaisala or third parties. Use of the software is governed by license terms and conditions included in the applicable supply contract or, in the absence of separate license terms and conditions, by the General License Conditions of Vaisala Group.

This product may contain open source software (OSS) components. In the event this product contains OSS components, then such OSS is governed by the terms and conditions of the applicable OSS licenses, and you are bound by the terms and conditions of such licenses in connection with your use and distribution of the OSS in this product. Applicable OSS licenses are included in the product itself or provided to you on any other applicable media, depending on each individual product and the product items delivered to you.

Table of Contents

English.....	5
Deutsch.....	13
Français.....	23
Español.....	33
Português.....	43
中文.....	53
日本語.....	61

Setup Roadmap

Before installing software and hardware, refer to the *viewLinc User Guide* and your specific device documentation for important information about minimum system requirements, security certificates, and licensing additional features.

If you require assistance with the setup of your system, the Vaisala Technical Support team is ready to assist you, free of charge ¹⁾ at helpdesk@vaisala.com.

1) Software product purchase entitles you to basic support, which is limited to hardware and software installation assistance. Comprehensive support is available through the purchase of a support plan.

Setup Checklist

viewLinc Enterprise Server	
<input type="checkbox"/>	A system administrator is available to help me with network and server configuration.
<input type="checkbox"/>	A dedicated network server is allocated to run viewLinc Enterprise Server 24x7. I know its IP address and hostname, and it meets minimum installation requirements.
<input type="checkbox"/>	The server's network connection is to a secure network, and the network type has been appropriately set on the server (Private or Domain , not Public).
<input type="checkbox"/>	I know the address of the Network Time Protocol (NTP) server in the network.
<input type="checkbox"/>	My laptop or workstation meets minimum client requirements.
<input type="checkbox"/>	I know where the program files and viewLinc database should be stored. A backup system is in place to make sure the files are recoverable in case of server failure.
<input type="checkbox"/>	I have the required security certificate and key files, or have the information needed to generate a security certificate and key during installation.
<input type="checkbox"/>	I have the USB drive that contains the viewLinc setup files and license key.

Devices	
<input type="checkbox"/>	I have a plan that lists device serial numbers, mounting and installation information.
<input type="checkbox"/>	IP addresses are allocated for my HMT330 series transmitters and data loggers' network connections (DHCP address reservations or static IPs).
<input type="checkbox"/>	Planned locations of wired devices are within 180 cm (6 feet) distance of power outlets.
<input type="checkbox"/>	I know if drilling or door access is required for probes (for example, freezer installation).

Network	
<input type="checkbox"/>	Ethernet network drops with RJ45 jacks are installed and tested at each location where a wired network connection is needed.
<input type="checkbox"/>	A Power over Ethernet (PoE) injector or switch provides operating power at network drops where a PoE powered device will be installed.
<input type="checkbox"/>	(Optional) I have installed a supported modem for SMS notifications.
<input type="checkbox"/>	(Optional) I have an account set up with a voice/SMS web notification service provider (Twilio), and an internet-accessible port available.

Optional License Keys			
I have the license key(s) required to enable one or more of the following features:			
<input type="checkbox"/>	Voice web service	<input type="checkbox"/>	Vaisala OPC UA Server
<input type="checkbox"/>	SMS web service	<input type="checkbox"/>	Web API
<input type="checkbox"/>	Modbus devices (non-Vaisala)		

Install viewLinc Enterprise Server Software

The viewLinc Enterprise Server software is delivered on a USB drive. Your license keys are included with the drive.

This is a simple overview of a new viewLinc Enterprise Server installation. There are many different installation and upgrade scenarios, depending on the type of installation and existing viewLinc version. Refer to the complete installation and upgrade procedures in the *viewLinc User Guide* and be sure to take backups of any existing data.

- ▶ 1. Insert the USB drive into your computer. If you want to install viewLinc Enterprise Server on a remote server (not the local computer), copy the **viewLincInstall.exe** file from the USB drive to the target server.
2. Run **viewLincInstall.exe**. The installation wizard starts.
3. Follow the instructions on screen to complete the installation wizard. If an earlier version of the viewLinc software is already installed on the server, it will be upgraded to the new version, and some installation choices may be skipped or modified as necessary.

viewLinc Enterprise Server software is now installed. The software runs in the background when the server is on, even when the user interface is not open.

Install Devices

See device-specific documentation for detailed instructions.

If you are upgrading an existing system without adding new devices, you can skip this section.

Create a device installation plan to keep track of data logger serial numbers and installation locations.

Install Access Points

Set up your **VaiNet access points** so that:

- They are able to connect to the Ethernet network
- They have the IP address or hostname of at least one NTP server they can access
- Their VaiNet channel number is set
- They have the IP address or hostname of viewLinc Enterprise Server
- They are in installation mode so that they can connect new VaiNet data loggers to the system

If you have new **WLAN access points**, configure their network and wireless settings.

When your wireless access points are set up, install them at their planned locations.

Install Data Loggers and Transmitters

Set up your **data loggers** and **transmitters** so that they are ready for installation.

If you have **Vaisala vNet PoE Data Logger Interfaces** for your wired data loggers, follow vNet documentation to connect and configure them.

Mount the data loggers and transmitters at their planned locations.

You can install many other types of devices to a viewLinc Monitoring System. See *viewLinc User Guide* for instructions on how to make use of these devices:

- viewLinc device hosts for connecting devices on a remote location
- Measurement devices with Modbus RTU or Modbus TCP/IP output
- Signal towers for visual and audible alarm status recognition
- Remote display terminals for viewing data

Discover and Accept Devices

New devices (3)

Data loggers must be accepted or added to the system before their data can be collected and used.

- ▶ 1. Open the viewLinc Enterprise Server user interface by double-clicking the desktop icon , or by entering the server's IP address and port in a web browser's address field.
2. Log in to viewLinc Enterprise Server as the default administrator user (username, **admin**, password, **viewL!ncAdm1n**).

It is recommended that you change the default administrator password immediately after logging in to ensure site security. Remember to record the new password in a safe location.

3. If this is the first time logging in as the admin user, a brief setup tour starts automatically. Complete the tour to familiarize yourself with system setup.
4. When viewLinc detects new devices, **New Devices** text appears at the top of the screen. Click it to open the **New Devices** window.
5. Compare the list of new devices to your device installation plan to verify that all your devices have been found. In the **Accept** column, select each device you want to add.
6. Select **Save** to register the selected devices with your viewLinc Enterprise Server.
7. Open the **Sites Manager** window and select the **Hosts and Devices** tab. On the **Hosts and Devices** tree, verify that all of your devices are now listed.

- If none of your devices are listed, check that network traffic to ports used by viewLinc (for example, 12600, 771, and 6767) is not blocked by a firewall.
- If some devices are missing, make sure they are configured, powered on, and connected to the network. Refer to the device-specific documentation for troubleshooting assistance.
- If the automatic discovery is not finding all of the devices, refer to the *viewLinc User Guide* for instructions on how to add or discover devices manually.

Configure viewLinc Enterprise Server

After installing, accepting or adding devices, you are ready to configure the viewLinc Enterprise Server to meet site-specific monitoring needs.

Download the **Planning Worksheet** (pdf) from **Overview > Getting Started** to define:

- Names for your monitored sites (Zones and Locations)
- Users, their assigned groups, rights and permissions
- Alarm thresholds required for different Locations
- Alarm notification delivery methods (email/SMS/command/voice), and notification escalation paths

Configuration instructions for the following tasks are provided in the *viewLinc User Guide* and online Help, and several interactive  tours are available from the viewLinc desktop **Help >**

Tours menu:

- Create Zones and Locations ( Create a Zone/Create a Location)
- Add groups and users ( Create a Group/Add a User)
- Add alarm templates ( Add/Apply a Threshold Alarm Template; Add/Apply an Alarm Notification Template)
- Set system preferences (such as system language, measurement unit display, email/SMS/voice settings)
- Optional setup tasks:
 - Assign permissions to groups
 - Set up schedules for alarms and user notifications
 - Define views for important Zones or Locations
 - Create predefined comments
 - Add signal towers or remote displays
 - Set up voice and/or SMS web services²⁾
 - Add non-Vaisala Modbus devices²⁾
 - Add support for Vaisala OPC UA Server²⁾

For more information about configuring viewLinc, refer to the instructions in the *viewLinc User Guide*.

²⁾ Requires purchase of feature license.

Test the System


System Test: Complete a standard system test to ensure all components are installed and functioning as expected. Recommended for non-GxP installations. For testing instructions, see *viewLinc User Guide*.

Validation: For controlled GxP environments, Vaisala provides a set of comprehensive **Installation Qualification and Operational Qualification Protocol (IQOQ)** and **GxP qualification documents**, to ensure installation of viewLinc and related system components meets your alarm monitoring requirements. Contact Vaisala to purchase viewLinc validation documents.

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