

# Vaisala RFL100

## Wireless Data Logger for Continuous Monitoring Systems

Vaisala HUMICAP® Humidity and Temperature Logger		RFL100	1	2	3	4	5	6	7	8	9	10	11	12	13	14
1	Probes															
	<b>1 Probe</b>															
2	PROBE 1 - Type and calibration	Probe type	Calibration points not customer selectable													
	Temperature (4 mm)	TMP115 0.5 m	-90, -45, 0, +45, +90													
	Temperature (4 mm)	TMP115 3 m	-90, -45, 0, +45, +90													
	Temperature (4 mm)	TMP115 3 m	-196, -90, -45, 0, +45, +90													
3	<b>None</b>															
4	<b>None</b>															
5	<b>None</b>															
6	Sensor protection	<b>No Sensor Protection</b>														
		<b>For TMP115 wide range Temperature probe (4 mm)</b>														
7	<b>None</b>															
8	Region Radio Frequency	<i>North America and Canada (FCC + CSA)</i>														
		<i>Europe (ETSI)</i>														
9	Battery type	<b>Standard Alkaline 1,5 V AA size</b>														
		<i>Lithium 1,5 V AA size</i>														
10	Logger installation	<b>Mounting bracket</b>														
		<i>Mounting bracket with magnets</i>														
11	Probe Extension Cable	<b>No Extension Cable</b>														
		<b>For HMP115 probe - attached to logger</b>														
	<i>3 meter Extension Cable</i>	For HMP110 & TMP115 probes														
	<i>10 meter Extension Cable</i>	For HMP110 & TMP115 probes														
	<i>3 meter Flat Cable for installations via door seal</i>	For HMP110 & TMP115 probes														
	<i>3 meter Cable and 3 meter Flat Cable</i>	For HMP110 & TMP115 probes														
	<i>10 meter Cable and 3 meter Flat Cable</i>	For HMP110 & TMP115 probes														
12	Probe Mounting accessories	<b>No Mounting accessories</b>														
		<i>Metal Clamp brackets for fridge or freezer wall mounting for probe with extension cable</i>														
		<i>Flange mounting kit for probe feed through installation through a duct wall</i>														
13	Quick Guide	<b>No Guide</b>														
		<i>Multilingual Guide</i>														
TOTAL																
QTY																
TOTAL VALUE																

**Additional information:**

*For ACCREDITED Calibration Please Order Item: HELACCRCAL  
For SPECIAL Calibration Please Order Item: HELSPECIALCAL*

**Selections in bold are included in the prices of the basic versions.**  
*Selections in italic are available at an extra price.*

Example of order code with typical settings:

RFL100	1	A	0	N	0	A	0	A	1	A	0	N	0	N
--------	---	---	---	---	---	---	---	---	---	---	---	---	---	---