

Safety information summary prepared by Vaisala

More detailed information about the contents of the liquids, storage, first aid measures in case of exposure can be found in the Safety Data Sheets prepared by Cargille Laboratories and delivered with the package.

Refractive Index Liquids (RI Liquids) by Vaisala (Supplier: Cargille Laboratories)

RI Liquid	Hazard symbol	Hazard statements	Exposure controls/ personal protection	Other information
Series AAA (1.300-1.399)		H413; May cause long lasting harmful effects to aquatic life.	Avoid contact with the eyes and skin and do not inhale fumes. Wash hands after handling the liquid.	P273; Avoid release to the environment. P501; Dispose of contents/container in accordance with local regulations. Must not be disposed together with household waste. Do not allow product to reach sewage system.
Series AA (1.400-1.459)	GHS08	H304; May be fatal if swallowed and enters airways.	Avoid contact with the eyes and skin and do not inhale fumes. Wash hands after handling the liquid.	P501; Dispose of contents/container in accordance with local regulations. Must not be disposed together with household waste. Do not allow product to reach sewage system.
Series A (1.460-1.570)	GHS08 GHS09	H304; May be fatal if swallowed and enters airways H411; Toxic to aquatic life with long lasting effects.	Avoid contact with the eyes and skin and do not inhale fumes. Wash hands after handling the liquid.	P273; Avoid release to the environment. P501; Dispose of contents/container in accordance with local regulations. Must not be disposed together with household waste. Do not allow product to reach sewage system.

Related material is available on our website at