



Order form
Purchaser

Valid from Mar 1st 2014

V 1.1

1200 Temperature Logger

			1	2	3	4	5	6	7	8	
Vaisala Veriteq Data Logger Package			DL1200					N	N	N	Price
Validation Type	<i>Validation</i>	<i>VL</i>	V								
	<i>Non Validation</i>	<i>SP</i>	S								
Channel Function 1	<i>Internal Temp Probe (Standard L Range Only)</i>		1								
Length Ext 1	<i>Temperature S/S</i>	<i>25 Feet</i>	2								
	<i>Temperature Water Proof</i>	<i>25 Feet</i>	3								
Channel Function 2	Not used			N							
Length Ext 2	<i>Temperature S/S</i>	<i>25 Feet</i>		A							
	<i>Temperature Water Proof</i>	<i>25 Feet</i>		B							
Additional Calibration	None	None			N						
Points	<i>1 Additional Point</i>	<i>Specify Additional Point</i>	1								
	<i>Additional 1 week delivery</i>	<i>2 Additional Points</i>	2								
Calibration Interval	Recommended	12 Months			N						
	Customer Specified	3 - 120 Months			A						
Reserved 4	None						N				
Reserved 1	None							N			
Reserved 2	None								N		

Example of order code with typical settings

25 Feet = 7.6 Meters

DL1200 V 1 A N N N N N