

Indigo™ 201 Analog Output Transmitter

		1 2 3 4 5				
Vaisala Indigo™ 201 Analog Output Transmitter		INDIGO201		0	A	0
1	Display	None <i>With metric units</i> <i>With non-metric units</i>	0 1 2			
2	Analog outputs type	4...20 mA 0...20 mA 0...1 V 0...5 V 0...10 V 1...5 V		A B C D E F		
3	For future use	None		0		
4	For future use	None			A	
5	Customization	None				0

Selections in bold are included in the prices of the basic version.
Selections in italic are available at an extra price.

Examples of order codes with typical settings:

						1 2 3 4 5				
INDIGO201	1	A	A	A	0					

Accessories

Remote probe cable for Indigo, 1 m
Remote probe cable for Indigo, 3 m
Remote probe cable for Indigo, 5 m
Remote probe cable for Indigo, 10 m

Item code

INDIGOCABLE1M
INDIGOCABLE3M
INDIGOCABLE5M
INDIGOCABLE10M